

JACKSON REES SYSTEMS FOR PEDIATRIC

autoclavable

Jackson Rees System I

consisting of:

- Silicone Rebreathing Bag 500 ml with open tube connection REF 66-00-141
- Connector, straight 22 mm O.D. to 22 mm O.D. REF 60-12-110
- Silicone Tube REF 64-00-230
- T-Connector REF 67-20-005
- Swivel Connector, right angled REF 60-60-002
- Connector, straight 22 mm I.D. REF 60-30-002
- O₂-Silicone Tubing (Length: 180 cm) REF 68-070-20

REF 67-12-000



ditto,

but Silicone Rebreathing Bag 500 ml with 5 mm hole (REF 66-00-441)


REF 67-11-000

Jackson Rees System II

consisting of:

- Silicone Rebreathing Bag 500 ml REF 66-00-011
- Connector, straight 22 mm O.D. to 22 mm O.D. REF 60-12-110
- Silicone Tube REF 64-00-230
- Jackson Rees Valve REF 67-10-002
- Swivel Connector, right angled REF 60-60-002
- Connector, straight 22 mm I.D. REF 60-30-002
- O₂-Silicone Tubing (Length: 180 cm) REF 68-070-20

REF 67-12-100



Jackson Rees System III

consisting of:

- Silicone Rebreathing Bag 500 ml REF 66-00-011
- Connector, straight 22 mm O.D. to 22 mm O.D. REF 60-12-110
- Jackson Rees Valve REF 67-10-222
- Connector, straight 22 mm O.D. to 22 mm O.D. REF 60-40-005
- Silicone Tube REF 65-00-200
- T-Connector REF 67-20-005
- Swivel Connector, right angled REF 60-60-002
- Connector, straight 22 mm I.D. REF 60-30-002
- O₂-Silicone Tubing (Length: 180 cm) REF 68-070-20

REF 67-12-200



Jackson Rees System IV

consisting of:

- Silicone Rebreathing Bag 500 ml REF 66-00-011
- Connector, straight 22 mm O.D. to 22 mm O.D. REF 60-12-110
- Jackson Rees Valve REF 67-10-002
- Swivel Connector, right angled REF 60-60-002
- Connector, straight 22 mm I.D. REF 60-30-002
- O₂-Silicone Tubing (Length: 180 cm) REF 68-070-20

REF 67-12-300



BREATHING SYSTEMS FOR ADULTS

autoclavable

Waters / Mapleson C

consisting of:

- Silicone Rebreathing Bag 2300 ml
- Waters Valve with O₂-Connection
- Swivel Connector, right angled
- Connector, straight 22 mm I.D.
- O₂-Silicone Tubing (Length: 180 cm)

REF 66-00-013
REF 67-10-003
REF 60-60-002
REF 60-30-002
REF 68-070-20




REF 67-12-500

Magill / Mapleson A

consisting of:

- Silicone Rebreathing Bag 2300 ml
- Magill Valve without O₂-Connection
- Swivel Connector, right angled
- T-Connector 22 mm I.D.
- 16 mm Silicone Respiration Tubing (Length: 180 cm)

REF 66-00-013
REF 67-10-333
REF 60-60-002
REF 60-11-533
REF 69-16-180



REF 67-12-600


EXPIRATORY VALVES

autoclavable, Material: Aluminium

Waters Valve

REF 67-10-003


Box 1



Magill Valve

REF 67-10-333

Box 1




JACKSON REES EXPIRATORY VALVES

autoclavable, Material: Aluminium

REF 67-10-002

Box 1



REF 67-10-222

Box 1



AYRES T-PIECES

autoclavable, Material: PSU

REF 67-20-005

Box 1



REF 67-20-004

Box 1

