

CPAP-MASK

for single use

CPAP-Respiration Sets

consisting of CPAP-Mask, Fixation tape, PEEP-Valve and Adapter for 22 mm tube connection

Size	PEEP-Valve	REF	Box
medium	2.5-10 cmH ₂ O	35-45-003	1
medium	5-20 cmH ₂ O	35-45-001	1
large	2.5-10 cmH ₂ O	35-45-004	1
large	5-20 cmH ₂ O	35-45-002	1



CPAP-Mask

consisting of CPAP-Mask and Adapter (REF 60-50-007)

Size	REF	Box
medium	35-40-100	5
large	35-40-000	5



Adapter, straight
for single use

REF 60-50-007	Box 10
---------------	--------

29 I.D. / 32 O.D.



PEEP-Valve

for single use, Material: PC

PEEP-Valve	Color	REF	Box
2.5-10 cmH ₂ O	green	85-13-004	1
5-20 cmH ₂ O	blue	85-13-002	1

PEEP-Valve

autoclavable, Material: PSU

PEEP-Valve	Color	REF	Box
2.5-10 cmH ₂ O	green	80-13-004	1
5-20 cmH ₂ O	blue	80-13-002	1



GASTRO-LARYNGEAL TUBE (G-LT)

Extraglottic tube with ventilation- and endoscopy lumen for control of airway patency during gastrointestinal endoscopic procedures.

The G-LT is designed for obtaining and maintaining control of airway patency during medium to long-term complex gastrointestinal endoscopic procedures performed on adults under deep sedation or general anesthesia while maintaining spontaneous or assisted ventilation.

- Diagnostic ERC (Endoscopic Retrograde Cholangiography)
- Diagnostic and therapeutic ERCP (Endoscopic Retrograde Cholangiopancreatography) for pancreatic and biliary disease:
 - Brush cytology and biopsy
 - Endosonographic gastrointestinal diagnosis
 - Removal of bile and pancreatic duct stones
 - Papillotomy or duodenal ampulla dilation
 - Treatment of biliary strictures (endoprosthesis, stent placement)
 - Palliation of malignancy: obstructive jaundice (drainage decompression)
- Enteroscopy
- Percutaneous Endoscopic Gastrostomy (PEG), particularly in neurologic patients.

autoclavable, Material: Silicone

Patient	Size	Color code	REF	Box
Adult	≥ 155 cm	blue	32-90-004	1

