

## CUFF PRESSURE GAUGES

For inflation and pressure control on tracheal tubes with high volume, low pressure cuffs  
Complete with connecting tube (100 cm)

New: with vacuum valve for evacuation of low pressure cuffs (beginning of 2012)

**Monitor**  
Ø 68 mm scale

REF 54-05-000



**Universal**  
Ø 68 mm scale

- Scala with two green ranges
- for tracheal tubes (22-32 cmH<sub>2</sub>O)
  - for Laryngeal Tubes and Laryngeal Masks (60-70 cmH<sub>2</sub>O)

REF 54-07-000



**Pocket**  
Ø 50 mm scale

REF 54-04-000



**Variabel**  
Ø 50 mm scale

- with adjustable pointer
- Scale in cmH<sub>2</sub>O and mmHg

especially for silicone ET Tubes with higher gas permeability

REF 54-06-000



**Sensitive**  
Ø 50 mm scale

Version	REF
with hook	54-03-001
without hook	54-03-000



**Pediatric**  
Ø 50 mm scale

Pressure range of 0-60 cmH<sub>2</sub>O

REF 54-02-001



## CUFF CONTROLLER

### Automatic Cuff Pressure Gauge for tracheal tubes with high volume, low pressure cuffs

- allows to seal with lowest pressures to reduce the risk of mucosal ischemia and pressure necrosis
- set pressure is automatically maintained
- Compensates leaks (pressure drop > 1cmH<sub>2</sub>O) to reduce the risk of aspiration with the possibility of pneumonia
- Pressure increase due to the diffusion of anesthetic gases is released to the set value
- Alarm system in case of leaks or disconnection
- Auto-Set Function inflates cuff automatically to 25 cmH<sub>2</sub>O (adult) / 10 cmH<sub>2</sub>O (pediatric)
- Pressure changes through touch pads
- Large LCD Display is easy to read
- Battery back-up in case of power failure
- Pressure can be changed between cmH<sub>2</sub>O and mmHg
- Universal Clamp for fixation to the standard rail



**Cuff Controller**  
with battery, universal clamp and  
connecting tubing (200 cm)

Patient	REF
Adult pre-adjusted to 25 cm H <sub>2</sub> O	<b>55-13-500</b>
Child pre-adjusted to 10 cm H <sub>2</sub> O	<b>55-13-501</b>

## ACCESSORIES

### Connecting Tube

Length	REF	Box
100 cm	<b>54-05-112</b>	10
200 cm	<b>54-05-113</b>	10



### Connecting Tube with one-way valve

Length	REF	Box
200 cm	<b>54-05-114</b>	10

