

## FACE MASKS (SILICONE)

autoclavable

### Bi-Mask with silicone lip

Silicone mask with internal plastic shell (PSU)

Size	Patient	Color code	REF	Box
2	Child	green	35-65-222	1
3	Adult, small	yellow	35-65-223	1
4	Adult, medium	red	35-65-224	1
5	Adult, large	purple	35-65-225	1
6	Adult, extra large	blue	35-65-226	1



### Bi-Mask with inflatable cushion

Silicone mask with internal plastic shell (PSU)

Size	Patient	Color code	REF	Box
3	Adult, small	yellow	35-75-223	1
4	Adult, medium	red	35-75-224	1
5/6	Adult, large	purple	35-75-225	1



### Mask with plastic dome and inflatable cushion

Size	Patient	Color code	REF	Box
2	Child	green	35-70-222	1
3	Adult, small	yellow	35-70-223	1
4	Adult, medium	red	35-70-224	1
5	Adult, large	purple	35-70-225	1
6	Adult, extra large	blue	35-70-226	1



### Mask with plastic dome and silicone lip

Size	Patient	Color code	REF	Box
3	Adult, small	yellow	80-11-003	1
4	Adult, medium	red	80-11-004	1
5	Adult, large	purple	80-11-005	1



### Rendell Baker Masks

Size	Patient	Connection	REF	Box
0	Newborn	22 mm I.D.	<b>35-60-000</b>	10
1	Infant	22 mm I.D.	<b>35-60-001</b>	10
2	Child, small	22 mm I.D.	<b>35-60-002</b>	10
3	Child	22 mm I.D.	<b>35-60-003</b>	10



### Pediatric masks, round

Size	Patient	Connection	REF	Box
0	Newborn	15 mm O.D.	<b>80-11-000</b>	10
1	Infant	15 mm O.D.	<b>80-11-001</b>	10
2	Child	22 mm I.D.	<b>80-11-002</b>	10



## FACE MASKS (PVC)

for single use

### Mask with inflatable cushion

Size	Patient	Color code	Connection	REF	Box
1	Newborn	—	15 mm O.D.	<b>35-40-201</b>	10
2	Infant	—	15 mm O.D.	<b>35-40-202</b>	10
3	Child	yellow	22 mm I.D.	<b>35-40-203</b>	10
4	Adult, small	red	22 mm I.D.	<b>35-40-204</b>	10
5	Adult, medium	purple	22 mm I.D.	<b>35-40-205</b>	10
6	Adult, extra large	blue	22 mm I.D.	<b>35-40-206</b>	10



## ACCESSORIES

**Hook ring, Ø 33 mm**  
fits masks with plastic dome,  
autoclavable, Material: Metall



**REF 80-11-266** Box 10

**Hook ring, Ø 30-33 mm**  
fits Bi-Masks,  
autoclavable, Material: PSU



**REF 35-60-266** Box 10

**Fixation tape**  
autoclavable, Material: Silicone

Patient	REF	Box
Child	<b>35-70-155</b>	10
Adult	<b>35-70-255</b>	10



**Fixation tape**  
for single use, Material: Isoprene

Patient	REF	Box
Adult	<b>35-30-255</b>	10

