

MANUJET III

The Manujet III with Jet Ventilation Catheters is recommended for emergency rooms, crash carts, ambulances and operating rooms as it guarantees a quick and efficient oxygenation of a patient.

The Manujet III is supplied in a complete kit with injector and catheters. The simple mechanical system is lightweight, portable and immediately ready for use.

Manujet III

with 4 m pressure hose, Jet-Ventilations-catheters acc. to Ravussin (13 G, 14 G, 16 G), 100 cm connecting tube, Endojet Adapter (for ET Tube) und Bronchoscope Adapter

REF 30-01-003



JET-VENTILATION CATHETERS ACC. TO RAVUSSIN

Allow a safe and atraumatic transtracheal access for oxygenation/ventilation with a manual jet (Manujet III) or an automatic High Frequency Jet Ventilator.

for single use, sterile


Size	Emergency Set with syringe	Box	Anesthesia Set without syringe	Box
16 G	REF 30-02-918-1	1	REF 30-02-018-1	5
14 G	REF 30-02-914-1	1	REF 30-02-014-1	5
13 G	REF 30-02-913-1	1	REF 30-02-013-1	5




JET CATHETERS (ORAL), single lumen

for single use, sterile

Patient	Length	REF	Box
Baby and Child	25 cm	30-07-250-1	5
Adult	30 cm	30-07-300-1	5
Adult	40 cm	30-07-400-1	5




ACCESSORIES




Connecting Tube
Length: 100 cm

REF 54-05-111 Box 10



Bronchoscope Adapter
for rigid bronchoscope

REF 60-70-005 Box 10



Endojet Adapter
for connection to ET Tube,
length: 50 cm,
for single use, sterile

REF 30-03-111-1 Box 5