

LARYNGOBLOC®

Fiber-optic Laryngoscope for single use

New Concept: "One-piece" blade and handle (Monoblock-Design)

Material: Polypropylene

Blade with Handle:

- torsion- and break proof due to "one-piece-design"
- resistant, comparable to metal blades
- no risk of contamination
- gentle to teeth
- non-sensitive to cold temperatures

Light Block with LED light source:

- extreme bright light
- constant light quality
- energy saving and long life time



Size	Description	standard handle	short handle	Box
0	Oxford	REF 36-40-000	REF 36-45-000	20
0	Miller	REF 36-30-000	REF 36-35-000	20
1	Miller	REF 36-30-001	REF 36-35-001	20
2	Macintosh	REF 36-20-002	REF 36-25-002	20
3	Macintosh	REF 36-20-003	REF 36-25-003 (Slim Design*)	20
4	Macintosh	REF 36-20-004	REF 36-25-004 (Slim Design*)	20

*Slim Design

thinner blade, especially for difficult intubation (obesity, pregnancy, reduced mouth opening)

standard handle



short handle



Light Block
with LED light source (Battery-Type: LR6-1.5V-AA), reusable

Version	REF	Box
standard (with 2 batteries)	36-10-000	1
short (with 1 battery)	36-15-000	1



Emergency Set
consisting of Laryngobloc® with short handle
• 1 x Macintosh # 4
• 1 x Oxford # 0
• 1 x Light Block, short

REF 36-90-100	Box 1
----------------------	-------



Standard Set
consisting of Laryngobloc® with standard handle
• 1 x Macintosh # 2, 3, 4
• 1 x Miller # 1
• 1 x Light Block, standard

REF 36-90-200	Box 1
----------------------	-------

