

MANUAL SUCTION PUMP

An effective and powerful manual suction pump for the quick suction of fluids from the oro- and nasopharyngeal cavities in emergency situations.

- cost effective alternative to electrical and foot powered devices
- portable and lightweight
- always ready for use as independent from any energy source
- easy to squeeze with maximum performance
- fits easily in most bags
- disposable collection canisters for adult and children
- soft and atraumatic suction tips
- valve prevents backflow of liquids to the patient
- MRI proof
- 15 mm I.D. Anschluss
 - fits directly to the ET Tube
 - allows connection of included suction tips
 - allows connection of included catheter adapter for funnel shape suction catheters

Set

with canister for adult and child, 3x suction tips (large, medium, small) and catheter adapter (for funnel shape suction catheters)

| Color | REF |
|-------|-----------|
| blue | 30-55-500 |
| olive | 30-55-501 |



ACCESSORIES

Suction Catheter for single use

| Size | CH | REF | Box |
|--------|----|-----------|-----|
| small | 14 | 30-57-305 | 10 |
| medium | 28 | 30-57-609 | 10 |
| large | 36 | 30-57-812 | 10 |



Suction Cap for single use

| | |
|---------------|--------|
| REF 30-56-001 | Box 10 |
|---------------|--------|



Catheter Adapter for funnel shape suction catheters, for single use

| | |
|---------------|--------|
| REF 30-58-060 | Box 10 |
|---------------|--------|



Canister with Suction Cap for single use

| Patient | REF | Box |
|----------------|-----------|-----|
| Child (50 ml) | 30-56-050 | 10 |
| Adult (250 ml) | 30-56-250 | 10 |



Carry Bag fits pump, canisters and suction tips

| |
|---------------|
| REF 30-64-000 |
|---------------|

