

PRESSURE TRANSDUCER HOLDERS

with 3-Joint-System, Universal clamp and clamping piece for transducer plates

Pressure Transducer Holder

- available in two length

Length	REF
30 cm	70-131-07
50 cm	70-231-07



PRESSURE TRANSDUCER HOLDERS for Pressure Transducer Plates and Stop-Cocks

autoclavable

Stop-cock holder

REF 78-10-000



Standard holder

REF 78-20-000




Universal holder with one flexible joint

REF 78-30-000



Universal holder with two flexible joints

REF 78-40-000



PRESSURE TRANSDUCER PLATES

to connect to the pressure transducer holders
autoclavable



for 2 pressure transducers

REF 79-000-52



for 2 "BRAUN" pressure transducers

REF 79-000-72