


## 19 MM RESPIRATION TUBINGS FOR ADULTS

autoclavable, Material: Silicone


Length	Color: green	Color: transparent	Box
30 cm	REF 66-22-030	REF 66-22-030tr	1
40 cm	REF 66-22-040	REF 66-22-040tr	1
60 cm	REF 66-22-060	REF 66-22-060tr	1
80 cm	REF 66-22-080	REF 66-22-080tr	1
90 cm	REF 66-22-090	REF 66-22-090tr	1
110 cm	REF 66-22-110	REF 66-22-110tr	1
130 cm	REF 66-22-130	REF 66-22-130tr	1
150 cm	REF 66-22-150	REF 66-22-150tr	1
180 cm	REF 66-22-180	REF 66-22-180tr	1
200 cm	REF 66-22-200	REF 66-22-200tr	1



## 19 MM RESPIRATION TUBINGS FOR ADULTS, with angled cuff

autoclavable, Material: Silicone, Color: green

Length	REF	Box
25 cm	66-22-925	1
40 cm	66-22-940	1




## 16 MM "LIGHTWEIGHT TUBINGS" FOR ADULTS

autoclavable, Material: Silicone, Color: green

- light weight due to smaller tube and thinner coil
- high stability
- longer lasting than autoclavable plastic tubings

Length	REF	Box
30 cm	69-16-030	1
45 cm	69-16-045	1
60 cm	69-16-060	1
90 cm	69-16-090	1
110 cm	69-16-110	1
130 cm	69-16-130	1
150 cm	69-16-150	1
180 cm	69-16-180	1
200 cm	69-16-200	1



## ACCESSORIES

autoclavable

**Hook**  
for 9 mm / 19 mm  
Respiration Tubings  
Material: Silicone

REF 66-00-002 Box 5



**Connector, straight**  
to connect and extend  
Respiration Tubings with  
22 mm cuffs, Material: POM

REF 60-15-110 Box 1




## 9 MM RESPIRATION TUBINGS FOR PEDIATRICS

autoclavable, Material: Silicone, Color: green


### 10 mm Cuffs

Length	REF	Box
25 cm	66-10-025	1
45 cm	66-10-045	1
60 cm	66-10-060	1
85 cm	66-10-085	1
110 cm	66-10-110	1
130 cm	66-10-130	1
150 cm	66-10-150	1
200 cm	66-10-200	1




### 10 mm / 22 mm Cuffs

Length	REF	Box
25 cm	66-11-025	1
45 cm	66-11-045	1
60 cm	66-11-060	1
85 cm	66-11-085	1
110 cm	66-11-110	1
130 cm	66-11-130	1
150 cm	66-11-150	1
170 cm	66-11-170	1
200 cm	66-11-200	1



### 22 mm Cuffs

Length	REF	Box
25 cm	66-12-025	1
60 cm	66-12-060	1
85 cm	66-12-085	1
110 cm	66-12-110	1
130 cm	66-12-130	1
150 cm	66-12-150	1
170 cm	66-12-170	1
200 cm	66-12-200	1



## 9 MM "LIGHTWEIGHT TUBINGS" FOR PEDIATRICS

autoclavable, Material: Silicone, Color: green

- light weight due to smaller tube and thinner coil
- high stability
- longer lasting than autoclavable plastic tubings

Length	REF	Box
25 cm	69-10-025	1
45 cm	69-10-045	1
60 cm	69-10-060	1
90 cm	69-10-090	1
110 cm	69-10-110	1
130 cm	69-10-130	1
150 cm	69-10-150	1
170 cm	69-10-170	1
200 cm	69-10-200	1

