


SILICONE RESUSCITATION BAGS

with barcode and serial number
autoclavable

New: Intake Valve includes Reservoir Valve (beginning of 2012)



Patient	40 cmH ₂ O pressure relief	20 cmH ₂ O pressure relief	without pressure relief	with manometer	with pressure proved manometer
Adult 2000 ml	REF 82-10-100	—	REF 82-10-101	—	—
Child 500 ml	REF 80-10-200	—	REF 80-10-201	REF 80-10-251	—
Infant 250 ml	REF 80-10-300	REF 80-10-306	REF 80-10-301	REF 80-10-351	REF 80-10-352

ACCESSORIES

Reservoir Bag
for single use

Size	REF	Box
600 ml	80-12-020	5
1800 ml	80-12-015	5
2600 ml	80-12-010	5



PEEP-Adapter
autoclavable, Material: PSU

REF 80-13-003	Box 1
---------------	-------



PEEP-Valve
autoclavable, Material: PSU

PEEP-Valve	Connection	Color	REF	Box
2.5-10 cmH ₂ O	22 mm O.D.	green	80-13-001	1
2.5-10 cmH ₂ O	30 mm I.D.	green	80-13-004	1
5-20 cmH ₂ O	22 mm O.D.	blue	80-13-005	1
5-20 cmH ₂ O	30 mm I.D.	blue	80-13-002	1



22 O.D.

30 I.D.

SILICONE RESUSCITATOR SETS

supplied in a rigid transparent carrying case

Adult

- Resuscitation Bag 2000 ml
- Mask with inflatable cushion # 5
- Reservoir Bag 2600 ml
- O₂ -Tubing 200 cm

REF 82-10-190

Child

- Resuscitation Bag 500 ml
- Mask with inflatable cushion # 3
- Reservoir Bag 1800 ml
- O₂ -Tubing 200 cm

REF 82-10-290

Infant

- Resuscitation Bag 250 ml
- Pediatric Mask, round # 1
- Reservoir Bag 600 ml
- O₂ -Tubing 200 cm

REF 80-10-390



SILICONE RESPIRATION-CASE

For a quick use in emergency situations. Case includes basic equipment like resuscitator, masks and Guedel airways. Ideal for hospitals, paramedics, dentists, industry, etc.

Adult

- Resuscitation Bag 2000 ml
- Mask with inflatable cushion # 3, 4, 5
- Guedel Airway # 3, 4, 5

REF 30-92-100

Child

- Resuscitation Bag 500 ml
- Pediatric Mask, round # 1
- Masken with inflatable cushion # 2, 3
- Guedel Airway # 0, 1, 2

REF 30-92-200



Infant

- Resuscitation Bag 250 ml
- Pediatric Masks, round # 0, 1
- Mask with inflatable cushion # 2
- Guedel Airway # 00, 0, 1

REF 30-92-300

Wall Fixation
for Respiration Case

REF 30-91-001



PVC-RESUSCITATION BAGS

Non-Rebreathing Valve and Intake Valve can be disassembled for single use



Patient	40 cmH ₂ O pressure relief	without pressure relief
Adult 1800 ml	REF 85-10-100	REF 85-10-101
Child 500 ml	REF 85-10-200	REF 85-10-201
Infant 250 ml	REF 85-10-300	REF 85-10-301

COMPLETE SETS

supplied in a rigid transparent carrying case for single use

Adult

- Resuscitation Bag 1800 ml
- Mask with inflatable cushion # 5
- Reservoir Bag 2600 ml
- Reservoir Valve
- O₂-Tubing 200 cm

REF 85-10-190

Child

- Resuscitation Bag 500 ml
- Mask with inflatable cushion # 3
- Reservoir Bag 1800 ml
- Reservoir Valve
- O₂-Tubing 200 cm

REF 85-10-290

Infant

- Resuscitation Bag 250 ml
- Mask with inflatable cushion # 2
- Reservoir Bag 600 ml
- Reservoir Valve
- O₂-Tubing 200 cm

REF 85-10-390



ACCESSORIES

Reservoir Bag

for single use

Size	REF	Box
600 ml	80-12-020	5
1800 ml	80-12-015	5
2600 ml	80-12-010	5



Reservoir Valve

for single use



REF 85-12-030 Box 10

PEEP-Valve

for single use, Material: PC

PEEP-Valve	Connection	Color	REF	Box
2.5-10 cmH ₂ O	22 mm O.D.	green	85-13-001	1
2.5-10 cmH ₂ O	30 mm I.D.	green	85-13-004	1
5-20 cmH ₂ O	22 mm O.D.	blue	85-13-005	1
5-20 cmH ₂ O	30 mm I.D.	blue	85-13-002	1



Peep Adapter

for single use



REF 85-13-003 Box 1

PVC-RESUSCITATION BAGS SETS

Valves are glued to the Resuscitation Bag, therefore completely disposable single use

consisting of:

- PVC-Resuscitation Bag (40 cmH₂O pressure relief) with face mask
- Reservoir Bag
- O₂-Tubing 200 cm
- PEEP Adapter



Adult

REF 84-10-195 Box = 6 kits



Child

REF 84-10-295 Box = 6 kits



Infant

REF 84-10-395 Box = 6 kits

MOUTH-TO-MOUTH-RESUSCITATOR

for emergency during life saving mouth to mouth resuscitation for single use, Material: PVC

- no mouth to mouth contact during CPR
- Non-Rebreathing-Valve protects rescuer from contamination
- bacterial/viral filter protects from infections
- compact size



REF 80-61-000 Box 1