

TOURNIQUET

with dual channel with dual channel for Bloodless Field, I.V. Regional Anesthesia and Bilateral Surgery
for electricity 100-240 VAC, 50/60 Hz

TOURNIQUET 9000

Digital Display

- large color coded LED display is easy to read
- precise monitoring shows actual cuff pressure

Precision Pressure Regulator

- fast cuff inflation and deflation
- regulator is easier to handle than touch pads
- limits pressure to max. 600mmHg

Flush-Taster

- to check for bleeding
- to release drug slowly after I.V. Regional Anesthesia

Alarm System

- visual and audible alarm to detect leaks, pressure deviations and low battery voltage

Color coded Coil Extension Hoses

- with Positive Locking Connectors (PLC)
- easy attachment and detachment
- for safe and leak-free inflation

Automatic Timer

- automatic time recording after cuff inflation
- provides elapsed cuff inflation time

Stand with basket

- very stable
- easy patient transport
- roomy storage for cuffs



Table Unit with Carrying Handle, Stand and Basket

REF 13-22-900

TOURNIQUET CUFFS

Double cuffs for I.V.-Regional Anesthesia



Double Cuffs

Patient / Size	Color code	Fabric Cuff, washable	Slim-Cuff, washable	Dispo-Cuff, sterile, for single use
Child / 30 cm	blue	REF 20-50-710	REF 20-70-710	REF 20-30-710SLZ-1
Arm / 46 cm	red	REF 20-50-712	REF 20-70-712	REF 20-30-712SLZ-1
Leg / 61 cm	green	REF 20-50-722	REF 20-70-722	REF 20-30-722SLZ-1

Further Tourniquets and Cuffs are available on request.