

## ENDOFIX ORAL II

oral ET Tube fixation for Anesthesia and Intensive Care  
for single use



Patient	Color code	for ET Tube	with elastic necktape	with foam necktape	Box
Child	yellow	4.5 - 6.5 mm	<b>REF 40-9006-010</b>	<b>REF 40-9005-010</b>	10
Adult	white	6.5 - 8,5 mm	<b>REF 40-9006-030</b>	<b>REF 40-9005-030</b>	10
Adult, large	blue	8.5 - 10.0 mm	<b>REF 40-9006-040</b>	<b>REF 40-9005-040</b>	10
Adult, extra large	purple	10.0 - 11.0 mm	<b>REF 40-9006-050</b>	<b>REF 40-9005-050</b>	10

**Foam Cushion**  
for single use

Packaging	REF	Box
sterile	<b>40-8111-1</b>	30
clean	<b>40-8222</b>	40



## RESCUEFIX

oral Tube fixation for emergencies  
for single use

In emergency situations a quick and safe tube fixation is vital. The Rescuefix allows a fixation without the use of adhesives as the integrated clamp fixes the tube safely and quickly. The separable necktape is adjustable in length and therefore suitable for patients with stiff neck as well. The integrated biteblock protects from tube occlusion.



Patient	Color code	for ET Tube	REF	Box
Child	yellow	5.0 - 6.5 mm	<b>40-9005-111</b>	10
Adult	white	7.0 - 10.0 mm	<b>40-9005-333</b>	10



## ENDOFIX NASAL

nasal ET Tube fixation for Intensive Care  
for single use

Version	REF	Box
with foam necktape	<b>40-4005</b>	25
with elastic necktape	<b>40-4006</b>	25



## ENDOTAPE

oral / nasal ET Tube fixation  
for single use

Patient	Length	REF	Box
Child and Adult	54 cm (21")	40-8018-044	10



## TUBE HOLDER

oral / nasal ET Tube fixation  
for single use

Patient	Length	REF	Box
Child and Adult	54 cm (21")	40-8010-033	10
Adult	65 cm (26")	40-8010-044	10



## BABY ENDOFIX

oral / nasal ET Tube fixation for babies  
for single use

for ET Tube	Color code	REF	Box
2.0 mm	transparent	40-1005-020	25
2.5 mm	yellow	40-1005-025	25
3.0 mm	red	40-1005-030	25
3.5 mm	green	40-1005-035	25
4.0 mm	blue	40-1005-040	25



## BABY-CAP

for comfortable fixation of tubes and catheters  
for single use, Material: 100% Cotton

- with individual and variable fixation possibility (oral / nasal) for tubes and catheters
- head size adjustable additionally ( $\pm 2$  cm), also covers intermediate sizes
- simple and quick access to the head
- color code for identification of different sizes
- as option with Endofix for safe fixation of ET Tube



Size	Color code	Head Size	Cap, alone	Cap with Endofix	for ET Tube	Box
0	white	19 - 22 mm	REF 40-1900-000	REF 40-1905-000	2.0 mm	5
1	orange	23 - 28 mm	REF 40-1900-001	REF 40-1905-001	2.5 mm	5
2	yellow	28 - 32 mm	REF 40-1900-002	REF 40-1905-002	2.5 mm	5
3	red	33 - 38 mm	REF 40-1900-003	REF 40-1905-003	3.0 mm	5
4	green	38 - 43 mm	REF 40-1900-004	REF 40-1905-004	3.5 mm	5
5	blue	43 - 50 mm	REF 40-1900-005	REF 40-1905-005	4.0 mm	5

## GUEDEL AIRWAY

### oropharyngeal airway with ET Tube fixation possibility

- rail and lateral clamp design allow direct fixation of the ET Tube to the Guedel Airway
- hard biteblock but soft and atraumatic tube
- same lumen from proximal to distal end allows the use of largest possible suction catheters
- one side open with smooth inner surface facilitates suctioning and easy cleaning of the product

autoclavable, Material: PUR

Size	Color code	Length	Guedel Airway	Box
000	light blue	3.5 cm	<b>REF 40-40-007</b>	10
00	dark blue	5.0 cm	<b>REF 40-40-008</b>	10
0	black	6.0 cm	<b>REF 40-40-009</b>	10
1	white	7.0 cm	<b>REF 40-40-010</b>	10
2	green	8.5 cm	<b>REF 40-40-020</b>	10
3	yellow	9.0 cm	<b>REF 40-40-030</b>	10
4	red	10.0 cm	<b>REF 40-40-040</b>	10
5	purple	11.0 cm	<b>REF 40-40-050</b>	10



for single use, sterile  
Material: EVA

Size	Color code	Length	Guedel Airway	Guedel Airway with elastic necktape	Box
000	light blue	3.5 cm	<b>REF 40-20-007-1</b>	<b>REF 40-26-007-1</b>	30
00	dark blue	5.0 cm	<b>REF 40-20-008-1</b>	<b>REF 40-26-008-1</b>	30
0	black	6.0 cm	<b>REF 40-20-009-1</b>	<b>REF 40-26-009-1</b>	30
1	white	7.0 cm	<b>REF 40-20-010-1</b>	<b>REF 40-26-010-1</b>	30
2	green	8.5 cm	<b>REF 40-20-020-1</b>	<b>REF 40-26-020-1</b>	30
3	yellow	9.0 cm	<b>REF 40-20-030-1</b>	<b>REF 40-26-030-1</b>	30
4	red	10.0 cm	<b>REF 40-20-040-1</b>	<b>REF 40-26-040-1</b>	30
5	purple	11.0 cm	<b>REF 40-20-050-1</b>	<b>REF 40-26-050-1</b>	30

## BITEBLOCK

### oral Tube fixations for ET Tube, Laryngeal Tube and Laryngeal Mask

- tube clips into biteblock and is therefore protected from occlusion



autoclavable, Material: PUR

Size	fits for	Biteblock	Box
4	ET Tube 6.5-8.0 mm	<b>REF 40-70-040</b>	10
5	ET Tube 8.5-11.0 mm	<b>REF 40-70-050</b>	10
	Laryngeal Tube # 2 + 2.5		
	Laryngeal Mask # 2 + 2.5		
6	Laryngeal Tube # 3, 4, 5	<b>REF 40-70-060</b>	10
	Laryngeal Mask # 3, 4, 5		

for single use, sterile  
Material: EVA

Biteblock	Box	Biteblock with elastic necktape	Box
<b>REF 40-60-040-1</b>	25	<b>REF 40-66-040-1</b>	20
<b>REF 40-60-050-1</b>	25	<b>REF 40-66-050-1</b>	20
<b>REF 40-60-060-1</b>	25	<b>REF 40-66-060-1</b>	20

## TRACHEOSTOMY NECKTAPES

Tube fixations for tracheostomized patients  
for single use

### with velcro fastener

Material: velours coated Foam, Color: blue

Patient	Length	REF	Box
Infant / Child	355 mm (14")	40-0005-011	10
Child / Adult	470 mm (19")	40-0005-033	10
Adult, extra large	520 mm (21")	40-0005-044	10



### with velcro fastener

Material: cotton coated Foam, Color: skin colored

Patient	Length	REF	Box
Infant / Child	340 mm (13")	40-0008-011	10
Child / Adult	490 mm (20")	40-0008-033	10
Adult, extra large	590 mm (23")	40-0008-044	10



### separable with velcro fastener

Material: cotton coated Foam, Color: skin colored

Patient	Length	REF	Box
Child / Adult	400 mm (16")	40-0066-033	10



### separable with hook fastener

Material: cotton coated Foam, Color: skin colored

Patient	Length	REF	Box
Child / Adult	350 mm (14")	40-0099-033	10



## TRACHEOSTOMY DRESSING

for single use, sterile

- breathable, acts as a second skin
- absorbs moisture
- avoids pressure sores
- approved as wound dressing

Size: 6.5 x 9.0 cm

REF 603014	Box 10
------------	--------



## "EAR DRESS"

for use after ear operations

for single use, Material: cotton coated foam with velcro fastening, Farbe: skin colored

- avoids pressure sores
- comfortable to wear
- safely protects the wound area

Patient	REF	Box
Child / Adult	40-7000	5

