

Endoscopy Mask

Allow simultaneous ventilation and endoscopy of the patient to increase patient comfort and the working quality of the doctor.

for single use

Application:

- Fiberoptic Intubation
The ET Tube can be railroaded orally or nasally along the bronchoscope through the flexible membrane. After intubation the membrane is pushed inside which allows easy removal of the mask.
- Bronchoscopy
The membrane guarantees an ideal seal and therefore an efficient ventilation before, during and after bronchoscopy without time pressure.
- Gastroenterology
- Trans Esophageal Echocardiography

Size	Patient	Membrane	REF	Box
0	Newborn	2.0 mm	30-40-000	6
1	Infant	2.0 mm	30-40-111	6
3	Child	3.0 mm	30-40-333	6
3	Child	5.0 mm	30-40-335	6
5	Adult	5.0 mm	30-40-555	6
5	Adult	10.0 mm	30-40-777	6



Fixation Tape

for Endoscopy Mask

Harness

Material: silicone, latex free, reusable, autoclavable

Patient	REF	Box
Child	35-70-155	10
Adult	35-70-255	10



Hook Ring

made of PE, Ø 32 mm fits for Endoscopy masks

REF 30-40-266 Box 10



Fixation Tape

Material: foam

for Mask	REF	Box
# 1	35-40-622	10
# 3+5	35-40-633	10



Bronchoscope Airway with external biteblock

for single use, sterile

- protects expensive equipment from damage
- guides bronchoscope and ET Tube

Size	Length	REF	Box
0	36 mm	30-40-400-1	10
2	65 mm	30-40-420-1	10
4	85 mm	30-40-440-1	10

