

Jet Ventilation

Manujet III

The Manujet III with Jet Ventilation Catheters is recommended for emergency rooms, crash carts, ambulances and operating rooms as it guarantees a quick and efficient oxygenation of a patient.

With 4 m pressure hose, Jet Ventilation Catheters acc. to Ravussin (13 G, 14 G, 16 G), 100 cm connecting tube, Endojet Adapter (for ET Tube) and Bronchoscope Adapter

REF 30-01-003

The Manujet III is supplied in a complete kit with injector and catheters. The simple mechanical system is light-weight, portable and immediately ready for use



Connecting Tube

100 cm long with Luer Lock


REF 54-05-111 Box 100



Bronchoscope Adapter

for connection to the rigid bronchoscope, 15 mm I.D. swivelling, with Luer Lock

REF 60-70-005 Box 10



Endojet Adapter

for connection to the Endotracheal Tube, 15 mm I.D. swivelling, 22 mm O.D. with expiration holes, 50 cm long, adjustable catheter with Luer Lock, sterile packed

REF 30-03-111-1 Box 5



Jet Ventilation Catheter acc. to Ravussin

Allows a safe and atraumatic transtracheal access for oxygenation/ventilation with a manual jet (Manujet III) or an automatic High Frequency Jet Ventilator.

Anesthesia Set

sterile packed

Patient	Size	REF	Box
Infant	16 G	30-02-018-1	5
Child	14 G	30-02-014-1	5
Adult	13 G	30-02-013-1	5

Emergency Set

sterile packed, 1 piece per box and syringe

Patient	Size	REF	Box
Infant	16 G	30-02-918-1	1
Child	14 G	30-02-914-1	1
Adult	13 G	30-02-913-1	1

