

# Manual Suction Pump

An effective and powerful manual suction pump for the quick suction of fluids from the oro- and nasopharyngeal cavities in emergency situations.

## Complete Set

with canister for adult and child, 3x suction tips (large, medium, small) and catheter adapter (for funnel shape suction catheters). color: blue

REF 30-55-500

ditto but, color: olive (military)

REF 30-55-501

- no power required
- portable and lightweight
- easy to squeeze with the maximum performance
- fits easily in most bags
- disposable collection canisters for adults and children
- cost effective alternative to electrical and foot powered devices
- soft and atraumatic suction tips
- valve prevents backflow of liquids to the patient
- MRI-proof
- 15 mm I.D. Connection
  - fits directly on the ET Tube
  - allows connection of included suction tips
  - allows connection of included catheter adapter for funnel shape suction catheters



### 250 ml Canister

for adults, with cap

REF 30-56-250 Box 10

### Suction Cap

with 15 mm I.D. connection

REF 30-56-001 Box 10



### Catheter Adapter

for funnel shape suction catheters

REF 30-58-060 Box 10



### 50 ml Canister

for children, with cap

REF 30-56-050 Box 10

### Carry Bag

fits pump, canisters, suction tips

REF 30-64-000



### Suction Tip

large, 36 CH

REF 30-57-812 Box 10



medium, 28 CH

REF 30-57-609 Box 10



small, 14 CH

REF 30-57-305 Box 10



## Meconium Aspirator

Material: PC, for single use

- 15 mm I.D. connection allows direct connection to the ET Tube for tracheal suctioning of aspirated meconium
- funnel shaped connection allows connection of a standard suction catheter for suctioning of meconium in mouth and pharynx

REF 49-30-000

Box 10

15 mm I.D.  
+ connection for  
Suction Catheter



funnel shaped connection