

# Pelvic Sling and Tourniquets

## Pelvic Sling

The Pelvic Sling helps to fix and stabilize fractures of the pelvis at the accident scene

The cuff which is integrated in the pelvic sling exerts a lateral pressure to the pelvis. This helps to reduce life threatening inner bleeding and pain significantly.

**Standard**, for single use  
for pelvis circumference  
80-120 cm (32"-47")  
with hand inflator

REF 34-10-004

**Large**, for single use  
for pelvis circumference  
100-140 cm (40"-55")  
with hand inflator

REF 34-10-005

**Standard**, reusable  
for pelvis circumference  
90-130 cm (36"-51")  
with manometer and hand inflator

REF 34-10-003



**Manometer 0-200 mmHg**  
for Pelvic Sling

REF 34-18-201

## Tourniquet "Tourny I"

with stainless steel buckle

REF 20-30-000 Box 20



## Tourniquet "Tourny II"

with aluminum buckle

REF 20-32-000 Box 20



## Tourniquet "Tourny III"

with plastic buckle

REF 20-30-100 Box 10



## Tourniquet for Child

elastic tape with velcro fastening

REF 20-30-005 Box 10



## Tourniquet Strap

elastic strap 40 cm,  
25 pcs perforated on one roll

REF 20-30-400 Box 1 roll

