

# Resuscitation Bag, single use

## PVC Resuscitation Bag

for single use

### Adult 1800 ml

with 40 cmH<sub>2</sub>O pressure relief  
without pressure relief

REF	Box
85-10-100	1
85-10-101	1

### Child 500 ml

with 40 cmH<sub>2</sub>O pressure relief  
without pressure relief

REF	Box
85-10-200	1
85-10-201	1

### Infant 250 ml

with 40 cmH<sub>2</sub>O pressure relief  
without pressure relief

REF	Box
85-10-300	1
85-10-301	1



## Face Mask

with inflatable cushion, for single use

Size	Patient	Color	Connection	REF	Box
1	Neonate	-	15 mm O.D.	35-40-201	10
2	Infant	-	15 mm O.D.	35-40-202	10
3	Child	yellow	22 mm I.D.	35-40-203	10
4	Adult small	red	22 mm I.D.	35-40-204	10
5	Adult medium	purple	22 mm I.D.	35-40-205	10
6	Adult large	blue	22 mm I.D.	35-40-206	10



## PEEP Valve

for single use, MRI-proof

2.5-10 cmH<sub>2</sub>O, color: green

	REF	Box
22 mm O.D.	85-13-001	1
30 mm I.D.	85-13-004	1



5-20 cmH<sub>2</sub>O, color: blue

	REF	Box
30 mm I.D.	85-13-002	1
22 mm O.D.	85-13-005	1



## PEEP Adapter

for single use

REF 85-13-003 Box 1



## Reservoir Bag

for single use

Size	REF	Box
2600 ml	80-12-010	5
1800 ml	80-12-015	5
600 ml	80-12-020	5



## Reservoir Valve

for single use

REF 85-12-030 Box 10



## Complete Set

for single use, supplied in a rigid transparent carrying case

### Adult

- with 1800 ml Resuscitation Bag
- Air Cushion Face Mask # 5
- Reservoir Bag 2600 ml
- Reservoir Valve
- 2 m O<sub>2</sub>-Tubing

REF 85-10-190

### Child

- with 500 ml Resuscitation Bag
- Air Cushion Face Mask # 3
- Reservoir Bag 1800 ml
- Reservoir Valve
- 2 m O<sub>2</sub>-Tubing

REF 85-10-290

### Infant

- with 250 ml Resuscitation Bag
- Air Cushion Face Mask # 2
- Reservoir Bag 600 ml
- Reservoir Valve
- 2 m O<sub>2</sub>-Tubing

REF 85-10-390

