

Stylets, Introducers and Tube Exchangers

METTS (Muallem ET Tube Stylet)

for difficult intubation

single use, sterile

- core made of metal, malleable (maintains curvature)
- flexible pre-formed tip
- graduation marks

| Size | Length | for ET Tube | REF | Box |
|-------|--------|-------------|--------------------|-----|
| 8 CH | 40 cm | ≥ 3.5 mm | 33-08-400-1 | 5 |
| 12 CH | 65 cm | ≥ 5.0 mm | 33-12-650-1 | 5 |
| 14 CH | 65 cm | ≥ 6.0 mm | 33-14-650-1 | 5 |



METTI (Muallem ET Tube Introducer)

for difficult intubation and tube exchange

for single use, sterile

- core made of plastic, semi-rigid (flexible)
- flexible pre-formed tip
- graduation marks


| Size | Length | for ET Tube | REF | Box |
|-------|--------|-------------|--------------------|-----|
| 12 CH | 80 cm | ≥ 5.0 mm | 33-12-800-1 | 5 |
| 14 CH | 80 cm | ≥ 6.0 mm | 33-14-800-1 | 5 |



Muallem Pipe Stylet

metal tube to pre-form and guide METTI
Length: 36 cm, 5 mm I.D.

REF 33-30-014 Box 1



Introducer

for difficult intubation

for single use, sterile

- lumen-design, Universal Connector included
- rigid pre-formed tip
- graduation marks

| Size | Length | for ET Tube | REF | Box |
|-------|--------|-------------|--------------------|-----|
| 14 CH | 65 cm | ≥ 5.0 mm | 33-50-650-1 | 5 |



Intubation Catheter

for fiberoptically assisted intubation

for single use, sterile

- lumen-design, Universal Connector included
- graduation marks

| Size | Length | for ET Tube | for Fiberoptic Scope | REF | Box |
|-------|--------|-------------|----------------------|--------------------|-----|
| 19 CH | 56 cm | ≥ 7.0 mm | ≤ 4.0 mm | 33-60-560-1 | 5 |




Tube Exchanger

for extubation and exchange of Endotracheal Tubes

for single use, sterile

- lumen-design, Universal Connector included
- graduation marks

| Size | Length | for ET Tube | REF | Box |
|-------|--------|-------------|--------------------|-----|
| 11 CH | 80 cm | ≥ 4.0 mm | 33-45-800-1 | 5 |
| 14 CH | 80 cm | ≥ 5.0 mm | 33-46-800-1 | 5 |
| 19 CH | 80 cm | ≥ 7.0 mm | 33-47-800-1 | 5 |



Universal Connector

15 mm / Luer Lock for oxygenation,
fits all sizes with lumen-design
Material: ABS

REF 33-60-005 Box 10



Stylets

short, malleable intubation stylet

autoclavable up to 134°C

- soft, atraumatic tip

without connector

| Size | Length | for ET Tube | REF | Box |
|--------|--------|---------------|------------------|-----|
| small | 34 cm | 2.5 - 3.5 mm | 90-10-111 | 25 |
| medium | 38 cm | 4.0 - 4.5 mm | 90-10-222 | 25 |
| large | 45 cm | 5.0 - 6.5 mm | 90-10-333 | 25 |
| large | 45 cm | 7.0 - 11.0 mm | 90-10-444 | 25 |

with adjustable 15 mm I.D. connector

| Size | Length | for ET Tube | REF | Box |
|--------|--------|---------------|------------------|-----|
| small | 34 cm | 2.5 - 3.5 mm | 90-10-001 | 10 |
| medium | 38 cm | 4.0 - 4.5 mm | 90-10-002 | 10 |
| large | 45 cm | 5.0 - 6.5 mm | 90-10-003 | 10 |
| large | 45 cm | 7.0 - 11.0 mm | 90-10-004 | 10 |

