

## AUTOMATIC PRESSURE INFUSORS

For the use with pressure infusion cuffs for automatic irrigation purposes during minimal invasive surgery.

### Precision Pressure Regulator

- very precise pressure regulation
- maintains selected pressure and guarantees automatic compensation
- limits pressure to max. 300 mmHg

### Protection Cover

- protects the device from liquid and damage

### Coil Extension Hose

- firmly connected to the device
- with male Luer Lock connector

### Universal Clamp

- horizontal or vertical position for attachment to drip stand or rail

### Pressure Infusor 100

for compressed air  
with single channel, complete with  
pressure hose (400 cm / 13 ft.)

Version	REF
with universal clamp	52-13-100
with stand and basket	52-15-100



### Pressure Infusor 200

for compressed air  
with switch for the connection of two pressure infusion  
cuffs, complete with pressure hose (400 cm / 13 ft.)

Version	REF
with universal clamp	52-13-200
with stand and basket	52-15-200



### Pressure Infusor 200 ELC

for electricity 230 V  
with switch for the connection  
of two pressure infusion cuffs

Version	REF
with universal clamp	52-13-200 ELC
with stand and basket	52-15-200 ELC