

ROLL-ON CUFFS

allow fast and easy blood exsanguination within seconds (an alternative to Esmarch Bandage)
 autoclavable, Material: Silicone, Color: transparent

Roll-On Cuff Set, small complete in case

- consisting of:
- Roll-On Cuffs # 1-4
 - manometer with connecting tube
 - measuring tape
 - brake-block, large or small

REF 21-95-400

Roll-On Cuff Set, large complete in case

- consisting of:
- Roll-On Cuffs # 3-6
 - hand pump with manometer and connecting tube
 - measuring tape
 - brake-block, large or small

REF 21-95-401



Roll-On Cuffs

Size	Patient limb	Color code	REF	Box
0	16-21 cm	grey	21-95-710	1
1	20-28 cm	white	21-95-711	1
2	26-33 cm	light blue	21-95-712	1
3	31-40 cm	yellow	21-95-719	1
4	38-49 cm	red	21-95-722	1
5	47-60 cm	green	21-95-727	1
6	58-70 cm	dark blue	21-95-729	1



ACCESSORIES

Hand Pump
 with manometer and
 connecting tube
 (100 cm)

REF 21-95-350



Brake-Block
 reusable, Material: Silicone

Size	REF	Box
small	21-95-200	1
large	21-95-201	1



Manometer
 with connecting tube
 (100 cm)

REF 21-95-300



Measuring Tape
 color coded, for easy identification
 of patient extremity and cuff size


REF 21-95-111 Box 10



ESMARCH BANDAGES

reusable, Material: Polyisoprene (PI)

Width	Length	REF	Box
6 cm	5 m	20-20-060	5
8 cm	5 m	20-20-080	5
10 cm	5 m	20-20-100	5
12 cm	5 m	20-20-120	5
15 cm	5 m	20-20-150	5



for single use, Material: Polyisoprene (PI)

Width	Length	REF	Box
6 cm	4 m	20-21-060-1	10
8 cm	3 m	20-23-080-1	10
8 cm	4 m	20-21-080-1	10
10 cm	4 m	20-21-100-1	10
15 cm	4 m	20-21-150-1	10



THREE CHAMBER BLOOD EXSANGUINATION GARMENT

Due to the three chamber-system the garment is inflated automatically from distal to proximal. With the 3 position zipper the garment can be adjusted to fit different size extremities for efficient exsanguination of the patients limb. Sleeve and bladder are exchangeable.

Automatic exsanguination for use with VBM Tourniquets 2800, 5000 and 5800.

Arm

REF 20-13-341



Leg

REF 20-13-343

