

## VBM INFUSOR (500 ml / 1000 ml / 3000 ml)

Pressure Infusion Cuff for rapid pressure infusion and irrigation purposes.  
reusable

### completely transparent

- allows clear view on the fluid bag and the fluid level
- 360° visible

### Wrap-Around Design with hook

- fluid bag can be replaced quickly and easily
- consistent pressure on the complete fluid bag

### Manometer

- pressure proved
- precise pressure range from 0-300 mmHg

### Hand Inflator

- large volume allows quick inflation

### Three-Way Stop Cock

- for pressure control
- for quick deflation

### reusable

- easy to clean (PU material)
- latex free



## ORDER INFORMATION



Size	with hand inflator, stop cock and color coded pressure indicator	with hand inflator and manometer	with hand inflator and pressure proved manometer	Box
500 ml	REF 57-02-050	REF 57-06-050	REF 57-08-050	1
1000 ml	REF 57-02-100	REF 57-06-100	REF 57-08-100	1
3000 ml	REF 57-02-300	REF 57-06-300	REF 57-08-300	1

## DISPO INFUSOR (500 ml / 1000 ml / 3000 ml)

Pressure Infusion Cuff for rapid pressure infusion and irrigation purposes.  
for single use

### Transparent Netting

- allows clear view on the fluid bag and the fluid level

### Pocket Design with hanger

- fluid bag can be replaced quickly and easily

### Pressure Gauge

- color coded with overpressure valve
- 360° visible

### Hand Inflator

- large volume allows quick inflation

### Three-Way Stop Cock

- for pressure control
- for quick deflation

### for single use

- reduces the potential for cross-contamination
- latex free



## ORDER INFORMATION



Size	with female Luer Lock	with hand inflator, stop cock and color coded pressure indicator	with hand inflator and stop cock	Box
500 ml	REF 56-01-050	REF 56-02-050	REF 56-03-050	5
1000 ml	REF 56-01-100	REF 56-02-100	REF 56-03-100	5
3000 ml	REF 56-01-300	REF 56-02-300		1

## PRESSURE INFUSION CUFFS 500 ml

Pressure Infusion Cuff for rapid pressure infusion and irrigation purposes.  
reusable

### Transparent Netting

- allows clear view on the fluid bag and the fluid level

### Pocket Design with hanger

- fluid bag can be replaced quickly and easily

### Manometer

- pressure proved
- precise pressure range from 0-300 mmHg

### Hand Inflator

- large volume allows quick inflation
- with release valve for pressure control

### reusable

- bladder and sleeve are replaceable



## ORDER INFORMATION



with hand inflator

REF 51-05-050 Box 1



with hand inflator and  
manometer

REF 51-06-050 Box 1



with hand inflator and pressure  
proved manometer

REF 51-08-050 Box 1

## PRESSURE INFUSION CUFFS 1000 ml

Pressure Infusion Cuff for rapid pressure infusion and irrigation purposes.  
reusable

### Transparent Netting

- allows clear view on the fluid bag and the fluid level

### 3 Versions

- pocket design with zip ①
- pocket design with hook ②
- wrap-around design with hook ③

### Manometer

- pressure proved
- precise pressure range from 0-300 mmHg

### Hand Inflator

- large volume allows quick inflation
- with release valve for pressure control

### reusable

- bladder and sleeve are replaceable



## ORDER INFORMATION



Version	with hand inflator	with hand inflator and manometer	with hand inflator and pressure proved manometer	Box
pocket design with zip	REF 52-05-100	REF 52-06-100	REF 52-08-100	1
pocket design with hook	REF 53-05-120	REF 53-06-120	REF 53-08-120	1
wrap-around design with hook	REF 53-05-100	REF 53-06-100	REF 53-08-100	1

## PRESSURE INFUSION CUFFS 1500 ml

Pressure Infusion Cuffs for irrigation purposes during minimal invasive surgery.  
reusable

### Transparent Netting

- allows clear view on the fluid bag and the fluid level

### 2 Versions

- pocket design with zip ①
- pocket design with hook ②

### Manometer


- pressure proved
- precise pressure range from 0-300 mmHg

### Hand Inflator

- large volume allows quick inflation
- with release valve for pressure control

### reusable

- bladder and sleeve are replaceable



## ORDER INFORMATION



Version	with female Luer Lock	with hand inflator and pressure proved manometer	with hand pump and pressure proved manometer	Box
pocket design with zip	REF 52-01-150	REF 52-08-150	REF 52-30-150	1
pocket design with hook	REF 53-01-150	REF 53-08-150	REF 53-30-150	1

## PRESSURE INFUSION CUFFS 3000 / 5000 ml

Pressure Infusion Cuffs for irrigation purposes during minimal invasive surgery.  
reusable

### Transparent Netting

- allows clear view on the fluid bag and the fluid level

### 2 Versions

- pocket design with zip ①
- pocket design with hook ②

### Manometer

- pressure proved
- precise pressure range from 0-300 mmHg

### Hand Inflator

- large volume allows quick inflation
- with release valve for pressure control

### reusable

- bladder and sleeve are replaceable



## ORDER INFORMATION



Version	with female Luer Lock	with hand inflator and pressure proved manometer	with hand pump and pressure proved manometer	Box
pocket design with zip	REF 52-01-300	REF 52-08-300	REF 52-30-300	1
pocket design with hook	REF 53-01-300	REF 53-08-300	REF 53-30-300	1



## MANOMETER 300 mmHg

pressure proved

REF 20-17-301



REF 20-16-300



with hand inflator  
and release valve

REF 20-18-300



with hand inflator  
and release valve

REF 51-09-666



## HAND INFLATOR

with release valve

REF 51-10-666



with three-way stop cock

REF 51-03-666



with color coded pressure gauge  
and three-way stop cock

REF 51-02-666



## PUMP

Hand pump  
with release valve

REF 51-30-000



Foot pump  
with silicone tubing (200 cm)  
and release valve

REF 51-20-000



## THREE-WAY STOP COCK

Material: metal

REF 51-02-000

