

SINGLE CUFFS



Patient / Size	Color code	Fabric Cuff • washable	Silicone / Velcro Cuff • autoclavable • easy to clean and dry quickly	Silicone Cuff • autoclavable • easy to clean and dry quickly	Slim Cuff • washable • low profile	Dispo Cuff • for single use • low profile • sterile packed • Box of 10 pcs.
Infant / 20 cm	white	REF 20-54-700	REF 20-64-700	—	REF 20-74-700	REF 20-34-700SLZ-1
Child / 30 cm	blue	REF 20-54-710	REF 20-64-710	REF 20-94-710	REF 20-74-710	REF 20-34-710SLZ-1
Arm / 35 cm	yellow	REF 20-54-711	REF 20-64-611	REF 20-94-711	REF 20-74-711	REF 20-34-711SLZ-1
Arm / 46 cm	red	REF 20-54-712	REF 20-64-612	REF 20-94-712	REF 20-74-712	REF 20-34-712SLZ-1
Lower Leg, Arm / 46 cm	red	REF 20-54-512 *c	REF 20-64-512 *c	—	REF 20-74-715 *c	REF 20-34-715SLZ-1 *c
Leg / 61 cm	green	REF 20-54-722 REF 20-54-522 *c	REF 20-64-522 *c	REF 20-94-722	REF 20-74-722 *c	REF 20-34-722SLZ-1 *c
Leg / 76 cm	dark blue	REF 20-54-727 REF 20-54-527 *c	REF 20-64-527 *c	REF 20-94-727	REF 20-74-727 *c	REF 20-34-727SLZ-1 *c
Leg / 86 cm	brown	REF 20-54-728 REF 20-54-528 *c	REF 20-64-528 *c	REF 20-94-728	REF 20-74-728 *c	REF 20-34-728SLZ-1 *c
Leg / 107 cm	dark red	REF 20-54-729	—	REF 20-94-729	REF 20-74-729 *c	REF 20-34-729SLZ-1 *c

*c = contour shaped

DOUBLE CUFFS for I.V. Regional Anesthesia



Patient / Size	Color code	Fabric Cuff • washable	Silicone / Velcro Cuff • autoclavable • easy to clean and dry quickly	Silicone Cuff • autoclavable • easy to clean and dry quickly	Slim Cuff • washable • low profile	Dispo Cuff • for single use • low profile • sterile packed • Box of 10 pcs.
Infant / 20 cm	white	REF 20-50-700	REF 20-60-700	—	—	—
Child / 30 cm	blue	REF 20-50-710	REF 20-60-710	REF 20-90-710	REF 20-70-710	REF 20-30-710SLZ-1
Arm / 35 cm	yellow	REF 20-50-711	REF 20-60-711	REF 20-90-711	—	—
Arm / 46 cm	red	REF 20-50-712	REF 20-60-712	REF 20-90-712	REF 20-70-712	REF 20-30-712SLZ-1
Leg / 61 cm	green	REF 20-50-722	REF 20-60-722	REF 20-90-722	REF 20-70-722	REF 20-30-722SLZ-1
Leg / 76 cm	dark blue	REF 20-50-727	REF 20-60-727	—	—	—
Leg / 86 cm	brown	REF 20-50-728	—	—	—	—
Leg / 107 cm	dark red	REF 20-50-729	—	—	—	—

SINGLE CUFFS with dual port



Patient / Size	Color code	Fabric Cuff • washable	Silicone / Velcro Cuff • autoclavable • easy to clean and dry quickly	Slim Cuff • washable • low profile	Dispo Cuff • for single use • low profile • sterile packed • Box of 10 pcs.
Infant / 20 cm	white	REF 20-54-700ZD	REF 20-64-700ZD	REF 20-76-700	REF 20-36-700SLZ-1
Child / 30 cm	blue	REF 20-54-710ZD	REF 20-64-710ZD	REF 20-76-710	REF 20-36-710SLZ-1
Arm / 35 cm	yellow	REF 20-54-711ZD	REF 20-64-611ZD	REF 20-76-711	REF 20-36-711SLZ-1
Arm / 46 cm	red	REF 20-54-712ZD	REF 20-64-612ZD	REF 20-76-712	REF 20-36-712SLZ-1
Lower Leg, Arm / 46 cm	red	REF 20-54-512ZD *c	REF 20-64-512ZD *c	REF 20-76-715 *c	REF 20-36-715SLZ-1 *c
Leg / 61 cm	green	REF 20-54-722ZD REF 20-54-522ZD *c	REF 20-64-522ZD *c	REF 20-76-722 *c	REF 20-36-722SLZ-1 *c
Leg / 76 cm	dark blue	REF 20-54-727ZD REF 20-54-527ZD *c	REF 20-64-527ZD *c	REF 20-76-727 *c	—
Leg / 86 cm	brown	REF 20-54-728ZD REF 20-54-528ZD *c	REF 20-64-528ZD *c	REF 20-76-728 *c	REF 20-36-728SLZ-1 *c
Leg / 107 cm	dark red	REF 20-54-729ZD	—	REF 20-76-729 *c	REF 20-36-729SLZ-1 *c

*c = contour shaped

DOUBLE CUFFS with dual port for I.V. Regional Anesthesia




Patient / Size	Color code	Fabric Cuff • washable	Slim Cuff • washable • low profile	Dispo Cuff • for single use • low profile • sterile packed • Box of 10 pcs.
Child / 30 cm	blue	REF 20-50-710ZD	REF 20-71-710	REF 20-31-710SLZ-1
Arm / 46 cm	red	REF 20-50-712ZD	REF 20-71-712	REF 20-31-712SLZ-1
Leg / 61 cm	green	REF 20-50-722ZD	REF 20-71-722	REF 20-31-722SLZ-1

TOURNIQUET DRAPES

to seal against fluid accumulation underneath the cuff
for single use, Material: TPE (latex free)


Size	clean packed	Box	sterile packed	Box
Arm	REF 24-00-100	40	REF 24-00-100-1	10
Leg	REF 24-00-200	40	REF 24-00-200-1	10
Leg, extra large	REF 24-00-300	40	REF 24-00-300-1	10



TOURNIQUET SLEEVES

for wrinkle-free padding
for single use


Patient	Size	Color code	clean packed	sterile packed	Box
Infant	≤ 20 cm	white	REF 23-00-700	REF 23-00-700-1	10
Child	≤ 30 cm	light blue	REF 23-00-710	REF 23-00-710-1	10
Arm	≤ 35 cm	yellow	REF 23-00-711	REF 23-00-711-1	10
Arm, long	≤ 46 cm	red	REF 23-00-712	REF 23-00-712-1	10
Leg	≤ 61 cm	green	REF 23-00-722	REF 23-00-722-1	10
Leg, long	≤ 76 cm	dark blue	REF 23-00-727	REF 23-00-727-1	10
Leg, extra long	≤ 86 cm	brown	REF 23-00-728	REF 23-00-728-1	10
Leg, super long	≤ 107 cm	dark red	REF 23-00-729	REF 23-00-729-1	10



CUFF COVER

keeps the cuff free from soiling and allows sealing the cuff with adhesive tape without damaging the cuff
for single use

Size	clean packed	sterile packed	Box
small	REF 20-11-111	REF 20-11-111-1	50
medium	REF 20-11-333	REF 20-11-333-1	50
large	REF 20-11-222	REF 20-11-222-1	50



MEASURING TAPE

color coded, for easy identification of patient extremity and cuff size

Version	REF	Box
reusable	20-95-111	10
for single use	20-95-222	50

