

SURGIMED HOSPITAL SUPPLIES

P.O. Box 12345 • Hatfield • 0028
Head Office - Tel: 012 430 2050 • Fax: 012 342 3064


E-mail: surgimed@netactive.co.za • Website: www.surgimed.co.za

SURGEON'S SCRUB BRUSHES

(ANCHOR - USA)

Code 2000 Super Surgeon's Scrub Brush


- Greater scrub-up comfort and efficiency and better soap retention is provided in this, the most efficient surgeon's scrub brush ever made.
- Special nylon bristles securely anchored by noncorrosive nickel-silver pins, (Soft Texture)
- Lightweight, patented hollow back block with grooved sides assures firm grip.
- Total tufts 112 (7 x 16)
- 15 Extra tufts to clean nails.



Code 1000 Standard Scrub Brush

The standard of the trade for thirty years. Identical to the Super Scrub except it does not have the row of short nail bristles. (Medium Texture)


- Total tufts 112 (7 x 16)



Code 1001 Cub Scrub Brush

Open-end, channel-back block with 126 securely anchored tufts of pure nylon bristles. (Firm Texture). Superior to any surgeon's brush in this price range.

- Total tufts 126 (7 x 18)



SURGIMED

INSTRUMENT CLEANING BRUSHES

Cannulated Instruments cleaning Brushes

Incorporating heavy duty nylon bristle ideal for Laparoscopic & small lumen instruments

71-00-200	3mm diameter	30cm long
71-00-210	3mm diameter	40cm long
71-00-220	3mm diameter	60cm long
71-00-230	4.5mm diameter	30cm long
71-00-240	4.5mm diameter	40cm long
71-00-250	4.5mm diameter	61cm long
71-00-260	6mm diameter	30cm long
71-00-270	6mm diameter	40cm long
71-00-280	6mm diameter	61cm long


Standard Magill cleaning Brushes

Soft nylon bristle ideal for rubber ET Tubes, Laryngeal masks, silicone tubing & catheter mounts


71-00-290	4.5mm diameter	25cm long (239 S)
71-00-300	8.0mm diameter	32cm long (239 M)
71-00-310	12.5mm diameter	40cm long (239 L)




D	OINST	Small wire brush with plastic handle	Bristle Section	Overall Length
			55 x 12mm	125mm



F	F7495	Large wire brush with wooden handle	Bristle Section	Overall Length
			100 x 25mm	240mm
F	F4302	Large wire brush with plastic handle	Bristle Section	Overall Length
			100 x 25mm	240mm



1 Code no: 71-00-050
Standard instrument cleaning brush, wooden handle with new blue polyester bristles which keeps brush from "flattening" and proiongs effectiveness
Overall length: 28cm Bristle section 11.43cm




2 Code no: 71-00-060
Tapered instrument cleaning brush, wooden handle with new blue polyester bristles.



3 Code no: 71-00-090
Small cleaning brush, with cranked nylon handle blue polyester bristles. Ideal for instruments jaws and small instruments such as Burrs
Overall length: 19cm Bristle section: 3.81cm

4 Code no: 71-00-350
Pull through extra fine cleaning brush incorporating goats hair bristle section. Ideal for fine suction tubes.



SURGIMED

SURGIMED